


# PCR Cabinet

## THE IDEAL DEVICE FOR DNA AMPLIFICATION

Work volume is protected by a UV beam.

Three widths are available:

- PCR 6 : 60 cm
- PCR 9 : 90 cm
- PCR 12 : 120 cm




1-2 Stainless Steel Brackets

3-4 Lateral Glass Panel

5- Back Glass Panel

6-7 Stainless Steel Brackets

8- Inclined Front Panel



### Structure

Made in laminated plywood. The internal surface is protected with brushed aluminum, efficient against UV lighting.

The work volume is delimited with Security Glass Stand. This type of glass is efficient and guarantees a full protection from UV beam. The inclined front glass enables manipulator to work comfortably.

### OPTIONS

- Removable work surface
- Customized Stand
- Front panel
- Hour counter

### Equipment

Germicid UV Lamp (wave length: 253,7 nm). This lighting is controlled by Timer 0-60 minutes.

